


SKIN ARCHITECT FILL & CARE®

Filler + Biostimulator
2 in 1

LAYERS
OF BEAUTY



Hypodermis



CROSS-3000®
technology




1023



VOLUME
26 mg/ml

DENSITY
23 mg/ml

SUBLIME
20 mg/ml



Low Cross-Linking ≈2%

- Robust cohesiveness
- Stable volumization, reducing the risk of migrations



3000 kDa Molecular weight

- Restores facial volume
- Provides structural integrity
- Reinforced durability



FILL & CARE



CROSS-3000® technology

- Consistent and systematic application of a rigorous purification process
- Reduces the BDDE agent to nearly insignificant levels (<0.1 ppm)
- Progressive-sequential sterilisation process in all production phases
- Eliminating any possible bacterial contamination and minimising endotoxin levels (<=0.5 EU/ml).

3000 kDa Molecular weight

- Similar to endogenous hyaluronic acid
- High biocompatibility
- Biostimulation of fibroblasts and adipocytes, and cell regeneration





**Optimal rheology
+ Viscoelasticity
+ Cohesiveness
+ Plasticity**

- High cohesiveness, monophasic and homogeneous
- Immediate final effect
- Natural results



Skin Architect Filler are a new generation of dermal fillers based on very high molecular weight hyaluronic acid, **specifically designed to protect the patient's health**, reducing adverse effects, physical sequelae and toxic waste.



High Safety Profile

- Low BDDE concentration (<0.1 ppm)*
- Marginal endotoxin levels (≤0.5 EU/ml)**
- Minimal cytotoxicity



**Low Cross-Linking
≈2%**

- Robust cohesiveness
- Stable volumization, no migrations
- Minimizes rejection and inflammation



Inflammation



Risk of oedema and granulomas



*BDDE Industry Standard Specification <2 ppm
**Industry standard specification <20 EU/ml

Fills

+ Improves overall appearance

- Distinct rheological properties designed to restore facial volume and structural integrity.
- Delivers immediate and noticeable wrinkle reduction through precise dermal support.
- Ensures sustained and deep hydration by effectively binding and retaining water within the dermal layers.
- Enhances skin aesthetics, promoting a luminous, radiant, and healthy appearance.

Biostimulates



+ Elastin and collagen

- Promotes long-term enhancement of skin through fibroblast stimulation and increased extracellular matrix protein synthesis: collagen and elastin upregulation.³
- Optimizes facial contour definition by inducing extracellular matrix remodeling and structural reorganization.
- Contributes to homeostasis by reducing erythema and melanin.^{4, 5, 6}

Regenerates

+ Regeneration of adipose tissue

- Restores aged and damaged skin by enhancing dermal density and improving structural firmness.^{7, 8}
- Facilitates fibroblast proliferation and migration, supporting extracellular matrix regeneration.⁹⁻¹⁰
- Activates adipose-derived stem cells (ADSCs), stimulating adipocyte proliferation and differentiation to restore lost volume and naturally restructure tissues.¹¹⁻¹²



Confidential, All rights reserved by Dermatriz Limited

**Hong Kong & Macau
Sole Distributor:**



**LESS,
BUT
MORE**



+852 9183-7388



info@dermatriz.com



Office 1501, 15/F, Island Centre, No.1 Great George Street, CWB, HK


Concentration
20 mg / ml

Elasticity (G')
80–200



SUBLIME

Smooths and defines



Injection depth

- dermis
- mucosa of the lips


Concentration
23 mg / ml

Elasticity (G')
190–350



DENSITY

Moderates and corrects



Injection depth

- dermis
- mucosa of the lips
- subcutaneous tissue

Concentration
26 mg / ml

Elasticity (G')
300–550



VOLUME

Provides structure and volume



Injection depth

- subcutaneous tissue
- suprapariosteal region

Composition	Monophasic cross-linked hyaluronic acid
Syringe	1 ml
Molecular weight	3000 kDa
Needle	27 G / 13 mm
Cross-linking agent	BDDE
Degree of cross-linking	~2%
Durability	Up to 12 months of duration

Clinical evidence



Case 1



1 Session

Skin Architect Filler Density (0.4 ml on lips)
Skin Architect Filler Volume (0.4 ml on lips)

Case 2



1 Session

Skin Architect Filler Sublime (5 ml)
Skin Architect Filler Volume (2.5 ml)

Case 3



1 Session

Skin Architect Filler Sublime (0.5 ml on lips)

Case 4



1 Session

Skin Architect Filler Sublime (0.6 ml on lips)

Case 5



1 Session

Skin Architect Filler Volume (1.5 ml on lips)
Skin Architect Filler Density (0.5 ml on lips)

Case 6



Middle third treatment

Skin Architect Volume (2 ml)
Skin Architect Density (1 ml)
Lower third treatment: Jawline definition
Skin Architect Volume (1 ml)

Case 7



Middle third treatment

2 ml · Skin Architect Volume (2 ml)
Nasolabial fold treatment
Skin Architect Sublime (1 ml)